

**TECHNICAL DATA SHEET**  
**Latex, Non-Sterile Medical Examination Gloves**  
**which meet ASTM D 3578-19**

**Intended Use**

Non-sterile, ambidextrous, single use, medical examination gloves are intended to prevent transmission of a wide variety of diseases to both patients and health care personnel. Specifically, gloves are intended to act as a barrier between the health care provider and patients to prevent contamination between patient and caregiver.

Characteristic		Testing Methodology	Acceptable Limit	
Dimension & Tolerance	Width (mm)	ASTM D 3767	Extra-Small (XS)	60 - 80
			Small (S)	70 - 90
			Unisize	75 - 95
			Medium (M)	85 - 105
			Large (L)	101 - 121
			X-Large (XL)	110 - 130
	XX-Large (XXL)		120 - 140	
Length (mm)		XS & S & Unisize: ≥ 220 Others: ≥ 230		
Thickness (mm)	Palm		≥ 0,08	
	Finger			
Water Tightness Test (Freedom from Holes)		ASTM D 5151	As prescribed by Inspection Level and AQL	
Physical Requirement	Before ageing	Tensile Strength	≥ 18MPa	
		Ultimate Elongation	≥ 650%	
	After ageing	Tensile Strength	≥ 14 MPa	
		Ultimate Elongation	≥ 500%	
Powder Residue	Powder Free Gloves	ASTM D 6124	≤ 2,0 mg/glove, Inspection Level N=5	
	Powdered Gloves		≤ 10,0 mg/dm <sup>2</sup> , Inspection Level N=2	
Latex Protein	Aqueous Extractable Protein	ASTM D 5712	≤ 200 µg/dm <sup>2</sup>	
	Antigenic Protein Content	ASTM D 6499	≤ 10 µg/dm <sup>2</sup>	

**Additional Information**

Storage Condition	Country of Origin
Store in a cool, dry place. Avoid excessive heat (40°C, 104°F). Shield open box from direct sun and fluorescent light.	Malaysia

**Appendix 1 List of Item Codes and Attributes**

## Appendix 1 (List of items)

*AMD Medicom Inc.*

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
1163	Medicom SafeTouch "Rejuvenate" Powder-Free Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Natural	Textured, Contain Lanolin & Vitamin E	Non-Sterile
1188	Medicom SafeBasics "Easy Fit" Powder-Free Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Natural	Textured	Non-Sterile
1189	Medicom SafeBasics "Easy Fit" Powdered Latex Medical Examination Gloves	200 Gloves/Box, 10 Boxes/Case, 2000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Natural	N/A	Non-Sterile
1124EX	Medicom SafeTouch "Connect" Powder-Free Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Natural	Chlorinated, Textured	Non-Sterile
1125EX	Medicom SafeTouch "Connect" Powdered Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Natural	Smooth	Non-Sterile

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
1189EX	Medicom SafeTouch "Easy Fit" Powdered Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Natural	N/A	Non-Sterile