

TECHNICAL DATA SHEET Latex Dental Dam

Intended Use

A rubber dam is composed of a thin sheet of latex with a hole in the center intended to isolate a tooth from fluids in the mouth during dental procedures.

Characteristic	Specification		
Material Composition	100% Natural rubber latex		
Description	Thin sheet in the shape of a rounded rectangle.		
Physical Properties (Before aging)	Tensile Strength	≥ 24 MPa	
	Elongation at Break	≥ 700%	
Physical Properties (After aging)	Tensile Strength	≥ 19 MPa	
	Elongation at Break	≥ 630%	
Powder Residue	≤ 2,0 mg/dam		
Latex Protein	Aqueous Extractable Protein	≤ 200 µg/dm2	
	Antigenic Protein Content	≤ 10 µg/dm²	
Dimension	Size	Width (mm)	Length (mm)
	5" x 5"	127	127
	6" x 6"	152	152
Thickness	Type	Thickness (mm)	
	Thin	0,14	
	Medium	0,18	
	Heavy	0,22	

Additional Information

Storage Condition	Country of Origin
Store in cool and dry condition. Avoid direct sunlight.	Malaysia

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
1255-HB	SafeTouch Natural Rubber Latex Dental Dams	52 Dental Dams/Box, 6 Boxes/Case, 312 Dental Dams/case	127 x 127 mm (5" x 5")	Blue	Heavy, Unflavored	Non-Sterile
1255-HG	SafeTouch Natural Rubber Latex Dental Dams	52 Dental Dams/Box, 6 Boxes/Case, 312 Dental Dams/case	127 x 127 mm (5" x 5")	Green	Heavy, Mint Flavor	Non-Sterile
1255-MB	SafeTouch Natural Rubber Latex Dental Dams	52 Dental Dams/Box, 6 Boxes/Case, 312 Dental Dams/case	127 x 127 mm (5" x 5")	Blue	Medium, Unflavored	Non-Sterile
1255-MG	SafeTouch Natural Rubber Latex Dental Dams	52 Dental Dams/Box, 6 Boxes/Case, 312 Dental Dams/case	127 x 127 mm (5" x 5")	Green	Medium, Mint Flavor	Non-Sterile
1255-TB	SafeTouch Natural Rubber Latex Dental Dams	52 Dental Dams/Box, 6 Boxes/Case, 312 Dental Dams/case	127 x 127 mm (5" x 5")	Blue	Thin, Unflavored	Non-Sterile
1255-TG	SafeTouch Natural Rubber Latex Dental Dams	52 Dental Dams/Box, 6 Boxes/Case, 312 Dental Dams/case	127 x 127 mm (5" x 5")	Green	Thin, Mint Flavor	Non-Sterile
1266-HB	SafeTouch Natural Rubber Latex Dental Dams	36 Dental Dams/Box, 6 Boxes/case, 216 Dental Dams/case	152 x 152 mm (6 x 6")	Blue	Heavy, Unflavored	Non-Sterile
1266-HG	SafeTouch Natural Rubber Latex Dental Dams	36 Dental Dams/Box, 6 Boxes/case, 216 Dental Dams/case	152 x 152 mm (6 x 6")	Green	Heavy, Mint Flavor	Non-Sterile
1266-MB	SafeTouch Natural Rubber Latex Dental Dams	36 Dental Dams/Box, 6 Boxes/case, 216 Dental Dams/case	152 x 152 mm (6 x 6")	Blue	Medium, Unflavored	Non-Sterile

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
1266-MG	SafeTouch Natural Rubber Latex Dental Dams	36 Dental Dams/Box, 6 Boxes/case, 216 Dental Dams/case	152 x 152 mm (6 x 6")	Green	Medium, Mint Flavor	Non-Sterile
1266-TB	SafeTouch Natural Rubber Latex Dental Dams	36 Dental Dams/Box, 6 Boxes/case, 216 Dental Dams/case	152 x 152 mm (6 x 6")	Blue	Thin, Unflavored	Non-Sterile
1266-TG	SafeTouch Natural Rubber Latex Dental Dams	36 Dental Dams/Box, 6 Boxes/case, 216 Dental Dams/case	152 x 152 mm (6 x 6")	Green	Thin, Mint Flavor	Non-Sterile