
TECHNICAL DATA SHEET
Bib

Intended Use

Dental bibs are made of paper intended to absorb saliva or moisture from the oral cavity during dental procedures.

Characteristic	Specification		
Material Composition	Tissue Paper Polyethylene Film		
Description	Folded piece of paper with a plastic backing. Embossed horizontal pattern on the body of the bib and vertical embossed patter on the vertical borders. (Refer to appendix 1 specific attribute)		
Dimension	**Unfolded Size**	**Unfolded Width (cm)**	**Unfolded Length (cm)**
	33 cm x 45 cm	33	45

Measurement Diagram



Visual



Additional Information

Storage Condition	Country of Origin
Store in a cool dry place	China

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color
8292	Medicom Safebasics Dryback Bibs	125 Bibs/Bag, 4 Bags/Case, 500 Bibs/Case	33 cm x 45 cm	Black

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color
8293-EX	Medicom Safebasics Dryback Bibs	50 Bibs/Bag,10 Bags/Case, 500 Bibs/Case	33 cm x 45 cm	Black